
Using CIVCO Products with FUJIFILM SonoSite, Inc. Systems

User Guide

FUJIFILM
Value from Innovation

SonoSite

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or 1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution | Federal (United States) law restricts this device to sale by or on the order of a physician.

SonoSite Edge, M-Turbo, S Series, SonoSite SII, X-Porte, MicroMaxx, NanoMaxx, TITAN, SonoSite, the SonoSite logo are registered and unregistered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions.

CIVCO is a registered trademark of CIVCO Medical Instruments Co., Inc.

DICOM is a registered trademark of the National Electrical Manufacturers Association.

All other trademarks are the property of their respective owners.

Patents: US 8,858,436; US 8,861,822; US 8,147,408; US 7,588,541; CN 103237499B; DE 602004027882D1; FR 1552792; GB 1552792, and JP 5972258.

Part Number: P01765-24

Publication Date: November 2017

Copyright © 2017 FUJIFILM SonoSite, Inc. All Rights reserved.

Using CIVCO Products with FUJIFILM SonoSite, Inc. Systems

Contents

Introduction	1
CIVCO reference guides required	2
Conventions, symbols, and terms	2
Customer assistance and comments	2
Needle guides	3
Preparing for the procedure	7
Verifying the needle path	8
Performing the needle procedure	12
Turning on guidelines	12
Cleaning and disinfecting	14
Temperature and humidity limits	14

Introduction

This document is used in conjunction with the CIVCO reference guides.

- ▶ This document provides information on setting up the ultrasound system and verifying the needle-guide path prior to use and is included with each starter kit.
- ▶ The CIVCO reference guides provide information on attaching the bracket, attaching the needle guide, cleaning and disinfecting procedures, sterilization procedures, and additional safety and labeling information and is included with each starter kit and replacement kits.

Please read both documents before using the bracket and needle guide on the ultrasound system. Both documents are for a reader familiar with ultrasound and vascular access techniques or invasive procedures; they do not provide training in sonography, vascular access, or clinical practices.

WARNING

Before using the system and needle guide, you must have training in interventional procedures as dictated by current relevant medical practices as well as in proper operation of the ultrasound system and transducer. During vascular access, the potential exists for serious complications including, without limitation, the following: pneumothorax, arterial puncture, guidewire misplacement, and risks normally associated with local or general anesthesia and recovery.

The needle guide is for use on the following ultrasound systems:

FC1 ultrasound system	SonoSite X-Porte ultrasound system
SonoSite Edge ultrasound system	SonoSite Edge II ultrasound system
S Series ultrasound system	SonoSite SII ultrasound system
M-Turbo ultrasound system	M-OB/GYN Office ultrasound system
MicroMaxx ultrasound system	TITAN high-resolution ultrasound system
SonoSite 180PLUS ultrasound system	NanoMaxx ultrasound system
iLook 25 personal imaging tool	

See the ultrasound system Help (SonoSite X-Porte only) or user guide for information about the ultrasound system, transducers, additional safety, and intended uses.

CIVCO reference guides required

Refer to the CIVCO reference guides for instructions on attaching, cleaning, and disinfecting the brackets and needle guides.

Conventions, symbols, and terms

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.

Symbols and terms used on the system are explained in the ultrasound system Help (SonoSite X-Porte only) or user guide.

Customer assistance and comments

Questions and comments are encouraged. FUJIFILM SonoSite is interested in your feedback regarding the system and this document. Please call FUJIFILM SonoSite at 1-888-482-9449 in the U.S. Outside the U.S., call the nearest FUJIFILM SonoSite representative.

For questions and comments regarding the CIVCO reference guides, contact CIVCO Medical Instruments at www.civcomedical.com, call (800) 445-6741 in the U.S., or call your local representative.

Needle guides

WARNINGS

- ▶ To avoid device damage or patient injury, do not use the P10, P10x, P10xp, P17, rP19x, rP19xp, P21x, or P21xp needle-guide bracket on patients with pacemakers or medical electronic implants. The needle-guide bracket for these transducers contains a magnet that is used to ensure the bracket is correctly oriented on the transducer. The magnetic field in direct proximity to the pacemaker or medical electronic implant may have an adverse effect.
- ▶ System software biopsy guidelines do not support the CIVCO Infiniti Needle Guidance system and should not be used together.

Needle guides assist in guiding a needle. The echoes of the anatomical target and the needle appear on the system display and assist in guiding the needle to the target.

When a needle guide-capable transducer is connected to the system, the system generates a pair of guidelines that represent the anticipated path of the needle. You can use the guidelines and depth cursor to assist in guiding the needle. See **“Turning on guidelines”** on page 12.

The guidelines display on the screen is intended for use with brackets that have installed needle guides only.

Caution

The guidelines display is not available for the SonoSite X-Porte ultrasound system when the centerline graphic is on. The centerline graphic aligns with the center mark of the transducer (if present) and serves as a reference mark for the center of the displayed image. The centerline graphic is not intended for use with needle-guide brackets or with any bracket-associated guidelines. For more information about the centerline graphic, see the *SonoSite X-Porte User Guide*.

The following table lists transducers that are needle guide-capable. Available needle guides may include the following types: fixed angle, multiple fixed angle, variable angle, and transverse angle.

Table 1: System and Transducer Compatibility

Transducer	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							

- 1 The starter kit for the L25 transducer series contains sterile needle guides in various gauges and lengths. See the CIVCO website at www.civcomedical.com for details.
- 2 The L25 transducer can be used with either 18- or 21-gauge needles and only supports a transverse needle guide (out-of-plane). Refer to the appropriate CIVCO reference guide for which needle length to use to ensure access to the desired target depth.

Table 1: System and Transducer Compatibility

Transducer	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						

1 The starter kit for the L25 transducer series contains sterile needle guides in various gauges and lengths. See the CIVCO website at www.civcomedical.com for details.

2 The L25 transducer can be used with either 18- or 21-gauge needles and only supports a transverse needle guide (out-of-plane). Refer to the appropriate CIVCO reference guide for which needle length to use to ensure access to the desired target depth.

Table 1: System and Transducer Compatibility

Transducer	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
L38xp	✓	✓											
P10									✓				
P10x			✓	✓	✓	✓	✓						
P10xp		✓											
P17									✓				
rP19x				✓		✓							
rP19xp		✓											
P21x			✓		✓		✓						
P21xp	✓	✓											

- 1 The starter kit for the L25 transducer series contains sterile needle guides in various gauges and lengths. See the CIVCO website at www.civcomedical.com for details.
- 2 The L25 transducer can be used with either 18- or 21-gauge needles and only supports a transverse needle guide (out-of-plane). Refer to the appropriate CIVCO reference guide for which needle length to use to ensure access to the desired target depth.

Preparing for the procedure

WARNINGS

To avoid or reduce the risk of infection to the patient:

- ▶ Use only FUJIFILM SonoSite or CIVCO approved needle guides, brackets, supplies, components, and accessories. Other brands may not properly fit FUJIFILM SonoSite transducers.
- ▶ Use only needle guides compatible with the systems and transducers listed in **Table 1, “System and Transducer Compatibility”** on page 4. Contact FUJIFILM SonoSite to determine compatible needle guides for a given system and transducer.
- ▶ Make sure that the bracket is attached correctly. Refer to the appropriate CIVCO reference guide.
- ▶ When using a multi-angle bracket, make sure that you select the same angle (A, B, or C) on both the bracket and the ultrasound system.
- ▶ Do not rely on the visibility of the needle tip on the system display. Use other tactile or visual indicators to determine whether you are at or in the vessel. (Example: indentation of anterior wall, decreased resistance as the needle enters the vessel lumen, or blood return in the needle.)
- ▶ Follow Universal Precautions when performing interventional procedures.

WARNING

Some transducer sheaths contain natural rubber latex and talc, which can cause allergic reactions in some individuals. Refer to 21 CFR 801.437, *User Labeling for Devices that Contain Natural Rubber*.

To prepare for the procedure

- 1 Verify that all components are not damaged.
- 2 Verify that you are using the correct bracket for the transducer.

The brackets for the P17, rP19 series, and P21 series transducers have a similar size and/or shape. (The bracket for the P21 series is labeled “P21”. The bracket for the rP19 series is labeled “P19”).

- 3 Confirm that the bracket is attached correctly.
 - ▶ Make sure that the bracket is secure. Apply pressure to snap it into place.
 - ▶ For the multi-angle bracket on the L38xi, L38xp, C60xf, C60xi, rC60xi, and C60xp transducers, make sure that the angle selected on the bracket matches the angle selected on the ultrasound system. See the ultrasound system Help (SonoSite X-Porte only) or user guide for information on system set up.
- 4 Verify the needle path. See **“Verifying the needle path”** on page 8.

5 In addition to the **Warnings** and **Cautions** in this document, observe the following precautionary measures:

- ▶ Needle guides must be attached over a sterile transducer sheath.
- ▶ FUJIFILM SonoSite does not recommend performing an anatomical survey for intracavity biopsies with the needle guide attached.
- ▶ Be aware that the needle may not be visible on the screen. Prior to the procedure, use the depth cursor to determine the target depth. Use this as guidance for depth of needle penetration.
- ▶ Use a straight, new needle and needle guide for each procedure.
- ▶ Refer to the appropriate CIVCO reference guide for information on attaching the bracket and needle guide.

Refer to the appropriate CIVCO reference guide for additional information to prepare for the procedure.

Verifying the needle path

Verify the needle path before the procedure to make sure that the relationship among the system on-screen guidelines, transducer, and needle guide is correct.

WARNING

To avoid injury to the patient, observe the following:

- ▶ Before performing needle path verification, check that all components are free of damage.
- ▶ Use a straight, new needle for each needle path verification procedure.
- ▶ Do not use the needle guide if the needle is not following the intended path.
- ▶ The needle used for this needle path verification must not be used for the actual procedure. Always use a new, sterile needle and needle guide for each procedure.

The following items are needed to verify the needle path (See **Figure 1, “L38 Transducer with the Bracket, Needle Guide, and Needle Attached”** on page 9):

- ▶ A needle guide-capable transducer
- ▶ Bracket (Non-disposable. The type of bracket you use depends on the transducer you are using.)
- ▶ Needle-guide kit (disposable)
- ▶ New, straight needle
- ▶ Ultrasound system

- ▶ Tank filled with test object fluid (All World Scientific code 3480) or distilled water. See [Figure 2, “L38 Transducer with Needle in Test Fluid”](#) on page 9


Figure 1 L38 Transducer with the Bracket, Needle Guide, and Needle Attached


Figure 2 L38 Transducer with Needle in Test Fluid

To verify the needle path

- 1 Check that all transducer, needle guide, and bracket components are free of damage.
- 2 Attach the bracket to the transducer, and attach the appropriate needle guide for the needle gauge. Refer to the appropriate CIVCO reference guide for instructions.
- 3 Set up the tank and transducer, and then insert the needle into the needle guide.

- 4 Immerse the transducer no more than 1/4-1/2 inch (6.35-12.7 mm) into a tank of test fluid or distilled water.
- 5 Connect the transducer, and turn on the ultrasound system.
- 6 Turn on the guidelines. See **“Turning on guidelines”** on page 12.
- 7 Move the needle down into the tank of test fluid or distilled water until its ultrasound image is visible on the system display (see **Figure 2, “L38 Transducer with Needle in Test Fluid”** on page 9).
- 8 Verify that the needle, as shown on the system display, is between the guidelines. See **Figure 3, “In-plane guides on the SonoSite X-Porte system”** on page 10 and **Figure 4, “In-plane guides on the MicroMaxx or TITAN system”** on page 11.

WARNINGS

- ▶ The guidelines are intended to provide only an indication of the expected path of the needle. Mechanical tolerances and needle flexibility are factors that may result in a needle diverging from the expected path.
- ▶ To avoid injury to the patient, if the needle enters from the unexpected side of the display, verify that the needle guide and bracket are correctly attached and that the orientation of the transducer is correct.
- ▶ If the needle is still not following the expected path between the guidelines, do not use the needle guide. Contact FUJIFILM SonoSite or your local representative.


Figure 3 In-plane guides on the SonoSite X-Porte system


Figure 4 In-plane guides on the MicroMaxx or TITAN system

Note:

The depth cursor measurement shown on the ultrasound system measures the distance from the tip of the needle to the origin of the needle guide. See [Figure 3](#) and [Figure 4](#).


Figure 5 Out of plane guide (all systems)

Performing the needle procedure

WARNING

To avoid injury or reduce risk of infection to the patient:

- ▶ Check that a needle path verification procedure has been performed before performing the needle procedure. For more information, see “[Verifying the needle path](#)” on page 8.
- ▶ Use a straight, new, sterile needle for each needle procedure.
- ▶ Sterilize reusable needle guides after each use.
- ▶ Thin needles can bend when entering tissue. Actual position must be verified by identifying the echoes from the needle.
- ▶ If the bracket, guide, or needle breaks during the procedure, verify that all pieces of equipment are recovered.
- ▶ If unable to visualize the needle at the target, ensure that the correct needle guide has been selected and that the bracket and needle guide are correctly attached.

Caution

To avoid damage to the transducer, FUJIFILM SonoSite recommends using Aquasonic gel. Using gels other than the one recommended by FUJIFILM SonoSite can damage the transducer and void the warranty. If you have questions about gel compatibility, contact FUJIFILM SonoSite or your local representative.

In addition to the **Warnings** and **Cautions** in this document, observe the following precautionary measures during a procedure:

- ▶ Use acoustic coupling gel for proper acoustic coupling. Apply a liberal amount of gel between the transducer and the body. Although most gels provide suitable acoustic coupling, some gels are incompatible with some transducer materials. A sample of Aquasonic[®] gel is provided in the needle-guide kit.
- ▶ Use sterile gel when performing a needle procedure.
- ▶ The guidelines are intended to provide only an indication of the expected path of the needle. Actual position must be verified by identifying the echoes from the needle.

Turning on guidelines

WARNING

Avoid using guidelines and depth indicators in **Zoom** or **Dual** mode. Viewing in **Zoom** or **Dual** mode can obscure important physical structures and reference points.

Guidelines are an optional feature to help with needle guidance. The availability of this feature depends on the type of transducer you are using. For transducers with a single-angle or multi-angle bracket, the touchpad moves the depth cursor.

To turn on the guidelines

- 1 On the ultrasound system, select the correct exam type.
- 2 Optimize the image.
- 3 In 2D imaging, do one of the following:
 - ▶ (SonoSite X-Porte) Tap **Needle Guide** and then tap **On**. If this control is hidden, tapping **More Controls** displays it.
Guidelines are not available in Zoom or Dual mode.
 - ▶ (M-Turbo, SonoSite Edge, SonoSite Edge II, S Series, and SonoSite SII) Select **Guide** on-screen.
 - ▶ (180PLUS) Press the **TEXT/PICTO** key. If pictograph is off, select the **Biopsy** icon to display the guidelines. If pictograph is on, select the **Forward** icon and select the **Biopsy** icon.
 - ▶ (FC1) On the touchscreen, tap **Needle Guide**.
 - ▶ (Other FUJIFILM SonoSite systems) Select **Biopsy** or **Guide** on-screen.
- 4 See the ultrasound system Help (FC1 and SonoSite X-Porte only) or user guide for more detailed instructions on basic operation and system settings.

To select the angle on a multi-angle bracket

- ❖ Tap Needle Guidelines and do one of the following:
 - ▶ (SonoSite X-Porte) Tap **Guide A**, **Guide B**, or **Guide C**.
 - ▶ (Other FUJIFILM SonoSite systems) Select **A**, **B**, or **C**.

To move the depth cursor

- ❖ Do one of the following:
 - ▶ (180PLUS) Use the trackball. If pictograph is on, press the  icon to return to the previous menu.
 - ▶ (SonoSite Edge, SonoSite Edge II) Press the buttons to increase or decrease depth.
 - ▶ (S Series, SonoSite SII) Turn the **Depth** knob.
 - ▶ (SonoSite X-Porte) On the control panel, use your finger to drag the depth cursor.

The target distance measurement value changes to reflect the distance between the needle guide origin and the depth cursor.

The origin is the point at which the needle enters the needle guide or needle insert. See [Figure 1](#).

To clear the guides

- ❖ Do one of the following:
 - ▶ (SonoSite X-Porte) Do either of the following tap **Needle Guide** and then tap **Off** (multi-angle), or
 - ▶ Tap **Needle Guide** (single-angle).
 - ▶ (M-Turbo, SonoSite Edge, SonoSite Edge II, S Series, and SonoSite SII) Do any of the following:
 - ▶ Select the angle again: **A**, **B**, or **C** (multi-angle).
 - ▶ Exit angle selection and press **Guide** (multi-angle).
 - ▶ Press **Guide** or **Biopsy** (single-angle).

Cleaning and disinfecting

Before cleaning, disinfecting, or sterilizing the bracket, refer to the appropriate CIVCO reference guide.

Temperature and humidity limits

For environmental limits, including shipping and storing, refer to the appropriate CIVCO reference guide.

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Schallkopf	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						

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Schallkopf	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xi				✓		✓							
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						

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Transductor	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						

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Transductor	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xi				✓		✓							
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						

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Sonde	FC1	SonoSite X-Porte	SonoSite Edge	SonoSiteEdgell	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						

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Sonde	FC1	SonoSite X-Porte	SonoSite Edge	SonoSiteEdgeII	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xi				✓		✓							
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						

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Trasduttore	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							

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Trasduttore	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						

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Transdutor	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							

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Transdutor	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						

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Transducer	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	SSeries	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						

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Transducer	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	SSeries	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xi				✓		✓							
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						

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Transduser	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							

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Transduser	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						

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Transduktor	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							

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Transduktor	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						

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Ηχοβολέας	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							

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Ηχοβολέας	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1,2}										✓	✓		✓
L25e ^{1,2}									✓				
L25x ^{1,2}			✓	✓	✓	✓	✓						
L25n ^{1,2}												✓	
L25xp ^{1,2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						
L38xp	✓	✓											

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Dönüştürücü	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							

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Dönüştürücü	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL38xp	✓	✓											
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1, 2}										✓	✓		✓
L25e ^{1, 2}									✓				
L25x ^{1, 2}			✓	✓	✓	✓	✓						
L25n ^{1, 2}												✓	
L25xp ^{1, 2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						

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轉換器	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
C8										✓	✓		
C8e			✓						✓				
C8x			✓		✓		✓		✓				
C35xp	✓	✓											
C60										✓	✓		
C60e									✓				
C60x			✓		✓		✓	✓					
C60xf	✓												
C60xi			✓		✓		✓						
C60xp		✓											
rC60xi				✓		✓							
HFL38									✓				
HFL38x			✓		✓		✓						
HFL38xi				✓		✓							
HFL38xp	✓	✓											

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轉換器	FC1	SonoSite X-Porte	SonoSite Edge	SonoSite Edge II	S Series	SonoSite SII	M-Turbo	M-OB/GYN Office	MicroMaxx	TITAN	SonoSite 180PLUS	NanoMaxx	iLook 25
HFL50x			✓	✓	✓	✓	✓						
HFL50xp	✓	✓											
ICT									✓	✓	✓		
ICTx			✓	✓	✓	✓	✓	✓					
ICTxp	✓	✓											
L25 ^{1,2}										✓	✓		✓
L25e ^{1,2}									✓				
L25x ^{1,2}			✓	✓	✓	✓	✓						
L25n ^{1,2}												✓	
L25xp ^{1,2}	✓	✓											
L38										✓	✓		
L38e									✓				
L38x			✓		✓		✓						
L38xi			✓	✓	✓	✓	✓						
L38xp	✓	✓											
P10									✓				
P10x			✓	✓	✓	✓	✓						

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