SonoSite Edge II

User Guide Errata

Manufacturer

FUJIFILM SonoSite, Inc.

21919 30th Drive SE

Bothell, WA 98021 USA

T: 1-888-482-9449 or 1-425-951-1200

F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.

Joop Geesinkweg 140

1114 AB Amsterdam,

The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd

114 Old Pittwater Road

BROOKVALE, NSW, 2100

Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, SonoSite Edge, and the SonoSite logo are registered and unregistered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions.

All other trademarks are the property of their respective owners.

Part Number: P19439-03

Publication Date: November 2017

Copyright © 2017 FUJIFILM SonoSite, Inc. All Rights reserved.

SonoSite Edge II Ultrasound System

User Guide Errata

Introduction	1
Document conventions	1
Getting Help	2
Saving images and clips	2
Measurement references	2
Measurement publications and terminology	2
Cleaning and disinfecting	5
Before getting started	6
Determining the required cleaning and disinfecting level	6
Clean and disinfect system and transducer to a high level (semi-critical uses)	7
Clean and disinfect system and transducer to a low level (non-critical uses)	13
Storing the transducer	16
Transporting the transducer	16

Introduction

Document conventions

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.
- ▶ Single-step procedures begin with ❖.

For a description of labeling symbols that appear on the product, see "Labeling Symbols" in the user guide.

Getting Help

For technical support, please contact FUJIFILM SonoSite as follows:

Phone (U.S. or Canada)	877-657-8118
Phone (outside U.S. or Canada)	425-951-1330, or call your local representative
Fax	425-951-6700
Email	ffss-service@fujifilm.com
Web	www.sonosite.com
Europe Service Center	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia Service Center	+65 6380-5581

Printed in the U.S.

The following content was either not included or incorrect in the system user guide; the revisions will be made in the next update.

Saving images and clips

The percentage icon in the system status area shows the percentage of space **available** in internal storage.

Measurement references

The following measurement references have been corrected.

Measurement publications and terminology

Effective Regurgitant Orifice (ERO) in mm²

Oh, J.K., J.B. Seward, and A.J. Tajik. *The Echo Manual*. 3rd Edition, Boston: Philadelphia: Lippincott Williams and Wilkins, (2007), 73-76, 210.

$$\text{ERO} = \text{MV Flow Rate} / \text{MR Vel} * 100$$

Ejection Fraction (EF), percent

Oh, J. K., J. B. Seward, A. J. Tajik. *The Echo Manual*. 3rd Edition, Philadelphia: Lippincott Williams and Wilkins, (2007), 115-116.

$$\text{EF} = ((\text{LVEDV} - \text{LVESV}) / \text{LVEDV}) * 100\%$$

where: EF = Ejection Fraction
LVEDV = Left Ventricular End Diastolic Volume
LVESV = Left Ventricular End Systolic Volume

Elapsed Time (ET) in msec

Oh, J. K., J. B. Seward, A. J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott Williams, and Wilkins, (2007), p.147, 148, figure 9-8

ET = time between velocity cursors in milliseconds

Left Atrial Volume

Lang R. et al. "Recommendations for Cardiac Chamber Quantification: A report from the American Society of Echocardiography's Guidelines and Standards Committee and the Chamber Quantification Writing Group, Developed in Conjunction with the European Association of Echocardiography, a Branch of the European Society of Cardiology." *J Am Soc Echocardiography* (2005), 18:1440-1463.

Lang R, Bierig M, Devereux R, et al. "Recommendations for Cardiac chamber Quantification by Echocardiography in Adults: An Update from the American Society of Echocardiography and the European Association of Cardiovascular Imaging." *J Am Soc Echocardiography* (2015), 28:1-39.

$$\text{LA Vol} = \pi/4(h) \Sigma(D1)(D2)$$

where: LA Vol = Left Atrial Volume in ml
h = Height of stacked oval disks making up the LA
D1 = Orthogonal minor axis
D2 = Orthogonal major axis

2-plane Simpson's rule (method of disks)

$$\text{LA Vol} = \pi/4(h) \Sigma(D1)(D2)$$

Simpson's algorithm divides the LA into a series of stacked oval disks where h is the height of the stacked disks and D1 and D2 are the orthogonal minor and major axes

1-plane Simpson's rule (method of disks)

$$LA\ Vol = \pi/4(h) \Sigma(D1)^2$$

Same as the 2-plane method of disks except there is an assumption that the stacked disks are circular.

The equation for the LA Vol Index is: LA Vol Index = LA Vol/BSA

MV Flow Rate in cc/sec

Oh, J.K., J.B. Seward, and A.J. Tajik. *The Echo Manual*. 3rd Edition, Philadelphia: Lippincott Williams and Wilkins, (2007), 73-76, 210.

$$Flow = PISA * Va$$

where: PISA = Proximal Isovelocity Surface Area

Va = aliasing velocity

Pressure Gradient (PGr) in mmHG

Oh, J.K., J.B. Seward, and A.J. Tajik. *The Echo Manual*. 3rd Edition, Philadelphia: Lippincott, Williams, and Wilkins, (2007), 63-66.

$$PGr = 4 * (Velocity)^2$$

Peak E Pressure Gradient (E PG)

$$E\ PG = 4 * PE^2$$

Peak A Pressure Gradient (A PG)

$$A\ PG = 4 * PA^2$$

Peak Pressure Gradient (PGmax)

$$PGmax = 4 * PV^2$$

Mean Pressure Gradient (PGmean)

Echocardiographic Assessment of Valve Stenosis: EAE/ASE Recommendations for Clinical Practice Journal of American Society of Echocardiography, (January, 2009), 4-5.

$$\text{PG mean} = \text{sum}(4v^2)/N$$

where: v = peak velocity at interval n
 N = the number of intervals in the Riemann sum.

Pressure Half Time (PHT) in msec

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), 73-74.

$$\text{PHT} = \text{DT} * 0.29$$

where: DT = deceleration time

Regurgitant Fraction (RF) in percent

Oh, J.K., J.B. Seward, and A.J. Tajik. *The Echo Manual*. 3rd Edition, Boston: Philadelphia: Lippincott Williams and Wilkins, (2007), 215-217.

$$\text{RF} = \text{RV} / \text{MV SV}$$

where: RV = Regurgitant Volume
 MV SV = Mitral Stroke Volume ($\text{Mitral CSA} * \text{Mitral VTI}$)

Mitral CSA = cross-sectional area calculated using Annulus Diameter

Regurgitant Volume (RV) in cc

Oh, J.K., J.B. Seward, and A.J. Tajik. *The Echo Manual*. 3rd Edition, Boston: Philadelphia: Lippincott Williams and Wilkins, (2007), 215-217.

$$\text{RV} = \text{ERO} * \text{MR VTI} / 100$$

Cleaning and disinfecting

This section updates the instructions for cleaning and disinfecting the ultrasound system and transducers.

Use the FUJIFILM SonoSite recommendations when cleaning or disinfecting your ultrasound system and transducer. Use the cleaning recommendations in the peripheral manufacturer's instructions when cleaning or disinfecting your peripherals.

See www.sonosite.com/transducers, for transducer images.

Notes

- ▶ The system and transducers must be cleaned and disinfected after each exam. It is important to follow these cleaning and disinfecting instructions without skipping any steps.
- ▶ For cleaning the TEEEx transducer, refer to the cleaning and disinfection instructions in the *TEEEx User Guide*.

Before getting started

- ▶ Follow the disinfectant manufacturer's recommendations regarding appropriate personal protective equipment (PPE), such as protective eyewear and gloves.
- ▶ Inspect the system and transducer to determine that they are free of any unacceptable deterioration, such as corrosion, discoloration, pitting, or cracked seals. If damage is evident, discontinue use, and contact FUJIFILM SonoSite or your local representative.
- ▶ Confirm that cleaning and disinfecting materials are appropriate for your facility's use. FUJIFILM SonoSite tests cleaners and disinfectants for use with the FUJIFILM SonoSite systems and transducers.
- ▶ Disinfectants and cleaning methods listed in this chapter are recommended by FUJIFILM SonoSite for efficacy and material compatibility with the products.
- ▶ Ensure that the disinfectant type and the solution strength and duration are appropriate for the equipment and application.
- ▶ Follow manufacturer recommendations and local regulations, when preparing, using and disposing of chemicals.

WARNINGS

- ▶ Ensure that cleaning and disinfecting solutions and wipes are not expired
- ▶ Some cleaners and disinfectants can cause an allergic reaction to some individuals

Caution

- ▶ Do not allow cleaning solution or disinfectant into the system connectors, or transducer connector.
- ▶ Do not use strong solvents such as thinner or benzene, or abrasive cleansers, since these will damage the exterior surfaces. Use only FUJIFILM SonoSite recommended cleaners or disinfectants.

Determining the required cleaning and disinfecting level

WARNING

The cleaning instructions contained in this chapter are based on requirements mandated by the American Food and Drug Administration (FDA). Failure to follow these instructions may result in cross contamination and patient infection.

The level of cleaning and disinfecting required for the system is dictated by the type of tissue it contacts during use. Use [Table 1](#) to determine the level of cleaning and disinfecting required.

Table 1: Choosing a cleaning and disinfecting method

Did any part of the system or transducer come in contact with broken skin, blood, mucosal membranes, or bodily fluids?	
<p>YES Came in contact with broken skin, blood, mucosal membranes, or bodily fluids.</p>	<p>Option A</p> <p>Go to “Clean and disinfect system and transducer to a high level (semi-critical uses)” on page 7.</p>
OR	
<p>NO Did not come in contact with broken skin, blood, or bodily fluids.</p>	<p>Option B</p> <p>Go to “Clean and disinfect system and transducer to a low level (non-critical uses)” on page 13.</p>

Spaulding classifications

Spaulding classifications (non-critical, semi-critical) determine the approach for cleaning and disinfecting medical equipment based on the device, the way it has been used, and the risk of infection. The system and transducers are designed for use within the Spaulding classifications of non-critical and semi-critical uses.

Option A Clean and disinfect system and transducer to a high level (semi-critical uses)

Use this procedure to clean and high-level disinfect the ultrasound system and transducer **whenever it has come into contact with blood, broken skin, mucosal membranes, or bodily fluids.**

Follow the manufacturer's instructions when using cleaners and disinfectants. The cleaners and disinfectants listed in the procedure are both chemically compatible and have been tested for efficacy with the system and transducers. Confirm that the cleaners and disinfectants are appropriate for your facility's use.

WARNINGS

- ▶ To avoid electrical shock, before cleaning, disconnect the system from the power supply.
- ▶ Wear the appropriate personal protective equipment (PPE) recommended by the chemical manufacturer, such as eyewear and gloves.

Caution

- ▶ Do not skip any steps or abbreviate the cleaning and disinfecting process in any way.
- ▶ Do not spray cleaners or disinfectants directly on the system surfaces or on system and transducer connectors. Doing so may cause solution to leak into the system, damaging it and voiding the warranty.
- ▶ Do not attempt to disinfect a transducer or transducer cable using a method or chemical that is not included here. This can damage the transducer and void the warranty.
- ▶ Use only FUJIFILM SonoSite recommended cleaners and disinfectants. Using a non-recommended disinfecting solution or incorrect solution strength can damage the System and transducer and void the warranty. Follow the disinfectant manufacturer's recommendations for solutions strengths.

Note

You must clean and disinfect both the ultrasound system and the transducer after each use, but only the transducer can be disinfected to a high level.

To clean and disinfect the system and transducer

- 1 Turn off** the system by pressing the **Power button**.
- 2 Unplug** the power cord from the outlet.
- 3 Remove** the disposable transducer sheath, if applicable.
- 4 Disconnect** the transducer from the system. Temporarily place it where it will not cross-contaminate clean equipment or surfaces while you clean the ultrasound console.
- 5 Clean** the exterior surfaces of the **ULTRASOUND SYSTEM** to remove any debris or bodily fluids. Use the following procedure:

- a Use either a pre-moistened wipe or a soft cloth dampened with cleaner or disinfectant. Choose a cleaner from the list of approved cleaners.

Approved cleaners/disinfectants for the ultrasound system

Cleaner/disinfectant	Minimum wet contact time ¹
SaniCloth AF3 ²	3 minutes
SaniCloth Plus	3 minutes
PI-Spray II	10 minutes

¹ For maximum effectiveness, the component being cleaned must remain wet with disinfectant for a minimum period of time.

² Qualified for use as an intermediate-level disinfectant for mycobacteria.

Refer to the cleaners and disinfection document available at www.sonosite.com for a more complete list of the most current cleaners and disinfectants.

- b Remove all gel, debris, and bodily fluids from the system.
 - c With a new wipe, clean the system, including the display, by wiping from clean areas to the soiled areas. This method helps to avoid cross-contamination.
 - d Observe the minimum wet contact time and refer to manufacturer's instructions. Monitor the system for wet appearance. Re-apply with a new wipe if no longer wet.
 - e Allow the ultrasound console to air dry in a clean, well-ventilated space.
- 6 Clean** the ultrasound system **STAND** to remove any debris or bodily fluids. For information about cleaning the stand, see the *Edge Stand and H-Universal Stand User Guide*
- 7 Clean** the **TRANSDUCER CABLE AND BODY** to remove any debris or bodily fluids. Use the following procedure:

- a Use either a pre-moistened wipe or a soft cloth dampened with cleaner or disinfectant. Choose a cleaner from the list of approved cleaners.

Approved cleaners/disinfectants for the transducer

Product	Compatible Transducer	Minimum wet contact time ¹
Sani-Cloth AF3 ²	C8x, C11x, C35x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, P10x, P11x ³ , rC60xi, rP19x	3 minutes
Sani-Cloth Plus	HFL50x, L52x ⁴	3 minutes
PI-Spray II	C8x, C11x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, L52x ⁴ , P10x, P11x ³ , rC60xi, rP19x	10 minutes

¹ For maximum effectiveness, the component being cleaned must remain wet with disinfectant for a minimum period of time.


² Qualified for use as an intermediate-level disinfectant for mycobacteria.

³ The L52x Transducer is for veterinary use only.

⁴The P11x Transducer is not licensed for use in Canada.

Refer to the cleaners and disinfection document available at www.sonosite.com for a more complete list of the most current cleaners and disinfectants.

- b Remove all gel, debris, and bodily fluids from the system.
- c With a new wipe, clean the cable and transducer, starting from the cable, wiping toward the scanhead. This method helps to avoid cross-contamination.

Caution | Do not allow moisture near the electronic components of the connector.

- d Observe the minimum wet contact time recommended by the manufacturer. Monitor the transducer for wet appearance. Re-apply with a new wipe if no longer wet.

8 Verify that all gel, debris, and bodily fluids have been removed from the system and transducer. If necessary, repeat steps 55, 6 & 7 with a new wipe.

WARNING Failure to remove all gel, debris, and bodily fluid could leave contaminants on the probe.

9 Examine the system, transducer and cable for damage, such as cracks or splitting where fluid can enter. If damage is evident, do not use. Instead, contact FUJIFILM SonoSite or your local representative.

10 Prepare the disinfectant for use.

- a Choose a high-level disinfectant from the list of approved disinfectants.

High-level disinfectants compatible with SonoSite Edge II transducers

Disinfectant	Compatible transducers	Temperature	Disinfectant Soak Duration
Cidex	C8x, C11x, C35x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, P10x, P11x ¹ , rC60xi, rP19x	25°C, 77°F	45 minutes
Cidex OPA	C8x, C11x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, L52x ² , P10x, rC60xi, rP19x	20°C, 68°F	12 minutes

¹ The P11x Transducer is not licensed for use in Canada.

² The L52x Transducer is for veterinary use only.

Refer to the cleaners and disinfection document available at www.sonosite.com for a complete list of the most current cleaners and disinfectants.

- b Check the expiration date on the bottle to ensure the disinfectant has not expired.
- c Mix or check that the disinfection chemicals have the concentration recommended by the manufacturer (for example, a chemical strip test).
- d Check that the temperature of the disinfectant is within the manufacturer's recommended limits.

11 Perform a high-level disinfection of the transducer. Use the following procedure:

- a Immerse the transducer in a high-level disinfectant solution.

WARNING High-level disinfectants can cause harm to the patient if not completely removed from the transducer. Follow the manufacturer's rinse instructions to remove chemical residue.

Caution

- ▶ Do not soak the transducer longer than recommended by the chemical manufacturer.
- ▶ Do not immerse the transducer connector in any disinfectant solution.
- ▶ Use only FUJIFILM SonoSite recommended cleaners and disinfectants. Using a non-recommended disinfecting solution or incorrect solution strength can damage or discolor the transducer and void the warranty.

- b** Ensure that the connector and most of the cable remain out of the fluid. You may submerge 2 inches (5 cm) of the cable proximal to the transducer scanhead.


Figure 1 Transducer component names

- 12 Rinse** the transducer **three separate times** using the following procedure:

WARNING

High-level disinfectants can cause harm to the patient if not completely removed from the transducer. Follow the manufacturer's rinse instructions to remove chemical residue.

- Rinse the transducer in clean, running water according to the disinfectant manufacturer's instructions (at least one minute).
- Ensure that the connector and at least 12-18 inches (31-46 cm) of the cable from the connector remain dry.
- Repeat this step until the transducer has been rinsed three separate times.

- 13 Dry** the transducer with a sterile lint-free cloth.

- 14 Dispose** of the disinfectant according to the manufacturer's guidelines.

15 Examine the transducer and cable for damage, such as cracks or splitting where fluid can enter.

If damage is evident, discontinue use of the transducer, and contact FUJIFILM SonoSite or your local representative.

Option **B** Clean and disinfect system and transducer to a low level (non-critical uses)

Use the following procedure to clean and disinfect the ultrasound system and transducer **if it has not come into contact with blood, broken skin, mucosal membranes, or bodily fluids.**

WARNING If the system or transducer has come into contact with any of the following, use the high-level cleaning and disinfection procedure. See **“Clean and disinfect system and transducer to a high level (semi-critical uses)”** on page 7:

- ▶ Blood
- ▶ Broken skin
- ▶ Mucosal membranes
- ▶ Bodily fluids

Follow the manufacturer's instructions when using cleaners and disinfectants. The cleaners and disinfectants listed in the procedure are both chemically compatible and have been tested for efficacy with the system and transducers. Confirm that the cleaners and disinfectants are appropriate for your facility's use.

WARNINGS

- ▶ To avoid electrical shock, before cleaning, disconnect the system from the power supply.
- ▶ Wear the appropriate personal protective equipment (PPE) recommended by the chemical manufacturer, such as eyewear and gloves.

Caution

- ▶ Do not spray cleaners or disinfectants directly on the system surfaces or on system and transducer connectors. Doing so may cause solution to leak into the system, damaging it and voiding the warranty.
- ▶ Use only FUJIFILM SonoSite-recommended cleaners and disinfectants. Using a non-recommended disinfecting solution or incorrect solution strength can damage the System and transducer and void the warranty. Follow the disinfectant manufacturer's recommendations for solutions strengths.
- ▶ Do not attempt to disinfect a transducer or transducer cable using a method or chemical that is not included here. Doing so can damage the transducer and void the warranty.

To clean and disinfect the system and transducer

- 1 **Turn off** the system by pressing the **Power button**.
- 2 **Unplug** the power cord from the outlet.
- 3 **Remove** the transducer sheath, if applicable.
- 4 **Disconnect** the transducer from the system. Temporarily place it where it will not cross-contaminate clean equipment or surfaces while you clean the ultrasound system.
- 5 **Clean** the exterior surfaces of the **ULTRASOUND SYSTEM** to remove any debris. Use the following procedure:
 - a Use either a pre-moistened wipe or a soft cloth dampened with cleaner and disinfectant. Choose a cleaner from the list of approved cleaners.

Approved cleaners/disinfectants for the ultrasound system:

Cleaner/disinfectant	Minimum wet contact time ¹
SaniCloth AF3 ²	3 minutes
SaniCloth Plus	3 minutes
PI-Spray II	10 minutes

¹ For maximum effectiveness, the component being cleaned must remain wet with disinfectant for a minimum period of time.

² Qualified for use as an intermediate-level disinfectant for mycobacteria.

Refer to the cleaners and disinfection document available at www.sonosite.com for a more complete list of the most current cleaners and disinfectants.

- b Remove all gel and debris from the system.
 - c With a new wipe, clean the system, including the display, by wiping from clean areas to the soiled areas. This method helps to avoid cross-contamination.
 - d Observe the minimum wet contact time and refer to manufacturer's instructions. Monitor the system for wet appearance. Re-apply with a new wipe if no longer wet.
- 6 **Clean** the ultrasound system **STAND** to remove any debris. For information about cleaning the stand, see the *Edge Stand and H-Universal Stand User Guide*.
- 7 **Clean** the **TRANSDUCER CABLE AND BODY** to remove any debris. Use the following procedure:

- a** Use either a pre-moistened wipe or a soft cloth dampened with cleaner or disinfectant. Choose a cleaner from the list of approved cleaners.

Approved cleaners/disinfectants for the transducer:

Product	Compatible Transducer	Minimum wet contact time ¹
Sani-Cloth AF3 ²	C8x, C11x, C35x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, P10x, P11x ³ , rC60xi, rP19x	3 minutes
Sani-Cloth Plus	HFL50x, L52x ⁴	3 minutes
PI-Spray II	C8x, C11x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, L52x ⁴ , P10x, P11x ³ , rC60xi, rP19x	10 minutes

¹For maximum effectiveness, the component being cleaned must remain wet with disinfectant for a minimum period of time.


²Qualified for use as an intermediate-level disinfectant for mycobacteria.

³The P11x Transducer is not licensed for use in Canada.

⁴The L52x Transducer is for veterinary use only.

Refer to the cleaners and disinfection document available at www.sonosite.com for a complete list of the most current cleaners and disinfectants.

- b** Remove all gel and debris from the transducer.
- c** With a new wipe, clean the cable and transducer, starting from the cable, wiping toward the scanhead. This method helps to avoid cross-contamination.

Caution | Do not allow moisture near the electronic components of the connector.

- d** Observe the minimum wet contact time recommended by the manufacturer. Monitor the transducer for wet appearance. Re-apply with a new wipe if no longer wet.
- 8 Verify** that all gel and debris have been removed from the system and transducer. If necessary, repeat steps 5, 6, & 7 with a new wipe.

9 Allow the ultrasound system, including the console and transducer, to air dry in a clean, well-ventilated space.

10 Examine the system, transducer and cable for damage, such as cracks or splitting where fluid can enter. If damage is evident, do not use. Instead, contact FUJIFILM SonoSite or your local representative.

Storing the transducer

To store the transducer

- 1 Make sure the transducer has been cleaned and disinfected as detailed in the previous section.
- 2 Store the transducer so that it hangs freely and vertically, and observe the following precautions:
 - ▶ Store the transducer away from any contaminated transducers.
 - ▶ Store the transducer in an environment that is safe and has good airflow. Do not store the transducer in closed containers or where condensation may occur.
 - ▶ Avoid direct sunlight and exposure to x-rays. Recommended storage temperature range is between 0° C (32° F) and +45° C (113° F).
 - ▶ If using a wall-mounted rack for storage, ensure that:
 - ▶ It is securely mounted.
 - ▶ The storage slots do not mar the transducer or cause damage to the cable.
 - ▶ The rack is sized and positioned to prevent the transducer from inadvertently falling.
 - ▶ Make sure the connector is supported and secure.

Transporting the transducer

When transporting the transducer, you must take precautions to protect the transducer from damage and avoid cross-contamination. Be sure to use a container approved by your organization.

To transport a soiled transducer for cleaning

A soiled transducer is one that has been contaminated and must be cleaned before using it in an exam.

- 1 Place the transducer in a clean, approved container.

WARNING To prevent cross-contamination or unprotected exposure of personnel to biological material, containers used to transport contaminated transducers should carry an ISO biohazard label similar to the following:

Caution Ensure the transducer is dry before placing it in a closed container. Condensation from a damp transducer can damage the connector.

- 2 Transport the transducer in the container to the point of processing. Do not open the container until the transducer is ready to be cleaned.

Caution Do not leave the transducer in a sealed container for long periods of time.

To transport a clean transducer

A clean transducer is one that has completed the cleaning and disinfection process, has been stored properly, and is ready to be used in an examination.

- 1 Place the transducer in a clean, approved container. To identify the transducer as clean, containers used to transport clean transducers should carry a cleanliness verification sticker or certificate.
- 2 Transport the transducer in the container to the point of use. Do not open the container until the transducer is ready to be used.

To ship a transducer

WARNING Whenever possible, avoid shipping a contaminated transducer. Before shipping, ensure the transducer has been cleaned and disinfected using the steps detailed in this chapter or according to special instructions received from FUJIFILM SonoSite. If you are returning the transducer to FUJIFILM SonoSite, document the disinfection on a "Declaration of Cleanliness," and attach it to the packing list.

- 1 Place the transducer in the shipping container and seal it. Do not allow any part of the transducer to protrude from the container.
- 2 Ship the transducer using the following precautions:
 - ▶ Clearly label the container as fragile.
 - ▶ Do not stack items on top of the shipping container.
 - ▶ Do not exceed the shipping temperature range: -35°C (-31°F) to $+65^{\circ}\text{C}$ (149°F).

- ▶ Do not open the shipping container until it reaches its final destination.
- ▶ After arrival, the transducer must be cleaned and disinfected before it can be used in an exam.

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Linkes Vorhofvolumen

Lang R. et al. „Recommendations for Cardiac Chamber Quantification: A report from the American Society of Echocardiography’s Guidelines and Standards Committee and the Chamber Quantification Writing Group, Developed in Conjunction with the European Association of Echocardiography, a Branch of the European Society of Cardiology.“ *J Am Soc Echocardiography* (2005), 18:1440-1463.

R. Lang, M. Bierig, R. Devereux et al. „Recommendations for Cardiac chamber Quantification by Echocardiography in Adults: An Update from the American Society of Echocardiography and the European Association of Cardiovascular Imaging.“ *J Am Soc Echocardiography* (2015), 28:1-39.

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Oh, J.K., J.B. Seward, and A.J. Tajik. *The Echo Manual*. 3rd Edition, Boston: Philadelphia: Lippincott Williams and Wilkins, (2007), 73–76, 210.

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Lang R. et al. «Recommendations for Cardiac Chamber Quantification: A report from the American Society of Echocardiography's Guidelines and Standards Committee and the Chamber Quantification Writing Group, Developed in Conjunction with the European Association of Echocardiography, a Branch of the European Society of Cardiology». *J Am Soc Echocardiography* (2005), 18:1440-1463.

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Lang R. et al."Recommendations for Cardiac Chamber Quantification: A report from the American Society of Echocardiography's Guidelines and Standards Committee and the Chamber Quantification Writing Group, Developed in Conjunction with the European Association of Echocardiography, a Branch of the European Society of Cardiology."*J Am Soc Echocardiography* (2005), 18:1440-1463.

Lang R, Bierig M, Devereux R, et al."Recommendations for Cardiac chamber Quantification by Echocardiography in Adults: An Update from the American Society of Echocardiography and the European Association of Cardiovascular Imaging."*J Am Soc Echocardiography* (2015), 28:1-39.

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Oh, J.K., J.B.Seward, and A.J.Tajik. *The Echo Manual*. 3rd Edition, Philadelphia: Lippincott Williams and Wilkins, (2007), 73-76, 210.

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Oh, J.K., J.B.Seward, and A.J.Tajik. *The Echo Manual*. 3rd Edition, Philadelphia: Lippincott, Williams, and Wilkins, (2007), 63-66.

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