**S SERIES ULTRASOUND SYSTEM**

*User Guide Supplement P15865-04*

---

**User Guide Supplement P15865-04**

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BROOKVALE, NSW, 2100

Australia

**Caution:** | United States federal law restricts this device to sale by or on the order of a physician.

S Series, V-Universal, SiteLink, SonoSite, and the SonoSite logo are trademarks or registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

DICOM is a registered trademark of the National Electrical Manufacturers Association.

P15865-04 06/2019

Copyright © 2019 FUJIFILM SonoSite, Inc. All Rights reserved.

# S Series Ultrasound System

## User Guide Supplement P15865-04

### Contents

Introduction .....	1
Intended uses .....	1
Saving images and clips .....	2
Connectivity setup .....	2
Maintenance .....	2
Troubleshooting .....	2
Electrical safety .....	3
Battery safety .....	3
Hazardous materials .....	3
Electromagnetic compatibility .....	3
Labeling symbols .....	12
Standards .....	13

This user guide supplement does the following:

- Provides additional safety information
- Updates the list of compatible accessories and peripherals

### Introduction

The user guide is intended for a reader familiar with ultrasound techniques and who has received training in sonography and clinical practices. Before using the S Series ultrasound system, you must receive such training.

### Intended uses

The clinician is positioned beside the patient and places the transducer onto (or into for invasive procedures) the patient's body where needed to obtain the desired ultrasound image.

## Saving images and clips

The percentage icon in the system status area shows the percentage of space used in internal storage. If you try to save an image or clip when no space remains, the system alerts you that internal storage is full. To resolve this issue, archive images and clips that you wish to keep, and then delete them from the system to free up space.

## Connectivity setup

For SiteLink™ issues, refer to the SiteLink Image Manager user guide. For DICOM® issues, such as storage commitment, archivers, and MPPS, refer to *Sending and Receiving DICOM Data on SonoSite Systems*.

## Maintenance

All maintenance requirements are described in this user guide. Performing maintenance procedures not described in the user guide may void the product warranty.

**WARNING:** | Do not modify the S Series ultrasound system.

## Troubleshooting

**Maintenance icon  appears on-screen.** Restart the system. If the issue recurs, system maintenance may be required. Note the number that appears in parentheses on the C: line, and contact SonoSite or your SonoSite representative.

**System prompts you to “ensure the USB device is valid.”** Make sure that the USB storage device does not have password protection enabled and is not defective.

Use the USB storage device included with the system.

**System prompts you to “ensure the USB device contains valid data.”** Make sure that the data is present on the USB storage device.

Reexport the original data onto the USB storage device.

Contact your system administrator.

**System displays the alert “Incompatible power supply. . .”** Use the power supply that shipped with the system. See [“Compatible accessories and peripherals”](#) on page 6.

**System displays the alert “Maximum number of procedure entries reached” when trying to create a patient information form.** Free up internal-storage space by archiving or exporting patient exams and then deleting them from the system.

**System displays the alert “Unable to save image or clip. You have reached the maximum number of images/clips allowed for a single patient.”** Delete any unwanted images or clips from the patient exam. For instructions, see the ultrasound system user guide.

## Electrical safety

- WARNING:** To avoid the risk of electrical shock:
- This equipment must be connected only to a supply mains with protective earth.
  - Do not let the bar code scanner touch the patient.
  - Do not let the power supply or connector touch the patient.
  - Do not touch the power supply and the patient at the same time.

## Battery safety

- WARNING:** Periodically, check to make sure that the battery charges fully. If the battery fails to charge fully, replace it.
- WARNING:** Do not ship a damaged battery without instructions from SonoSite Technical Support.

## Hazardous materials

- WARNING:** Products and accessories may contain hazardous materials. Ensure that products and accessories are disposed of in an environmentally responsible manner and meet federal and local regulations for disposing hazardous materials.

## Electromagnetic compatibility

- WARNING:** The S Series ultrasound system should not be used adjacent to or stacked with other equipment. If such use occurs, verify that the S Series ultrasound system operates normally in that configuration.

**Caution:** To avoid the risk of increased electromagnetic emissions or decreased immunity, use only accessories and peripherals recommended by SonoSite. Connection of accessories and peripherals not recommended by SonoSite to the ultrasound system may result in malfunction of the ultrasound system or other medical electrical devices in the area. Contact SonoSite or your local representative for a list of accessories and peripherals available from or recommended by SonoSite.

**Caution:** Electrostatic discharge (ESD), or static shock, is a naturally occurring phenomenon. ESD is common in conditions of low humidity, which can be caused by heating or air conditioning. ESD is a discharge of the electrical energy from a charged body to a lesser or non-charged body. The degree of discharge can be significant enough to cause damage to a transducer or an ultrasound system. The following precautions can help reduce ESD: anti-static spray on carpets, anti-static spray on linoleum, and anti-static mats.

## Separation distance

### Recommended separation distances between portable and mobile RF communications equipment and the S Series ultrasound system

The S Series ultrasound system is intended for use in an electromagnetic environment in which radiated radio frequency (RF) disturbances are controlled. The customer or the user of the S Series ultrasound system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the S Series ultrasound system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter Watts	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

*Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.*

## Compatible accessories and peripherals

SonoSite has tested the S Series ultrasound system with the following accessories and peripherals and has demonstrated compliance to the requirements of IEC60601-1-2:2007.

You may use these SonoSite accessories and third-party peripherals with the S Series ultrasound system.

**WARNING:** Use of the accessories with medical systems other than the S Series ultrasound system may result in increased emissions or decreased immunity of the medical system.

**WARNING:** Use of accessories other than those specified may result in increased emissions or decreased immunity of the ultrasound system.

### Accessories and peripherals compatible with S Series ultrasound system

Description	Maximum Cable Length
C11x transducer	6.5 ft/2.0 m
C60x transducer	6.0 ft/1.8 m
HFL38x transducer	6.0 ft/1.8 m
HFL50x transducer	6.0 ft/1.8 m
ICTx transducer	6.0 ft/1.8 m
L25x transducer	8.0 ft/2.4 m
L38xi transducer	5.5 ft/1.7 m
L52x transducer	7.9 ft/2.4 m
P10x transducer	6.5 ft/2.0 m
P21x Transducer	6.5 ft/2.0 m
SLAx transducer	8.0 ft/2.4 m
TEEx Transducer	7.5 ft/2.3 m
Bar code scanner	4.8 ft/1.5 m
Battery for PowerPack	—
Battery Pack	—
Battery PowerPack	—
Black & white printer	—

### Accessories and peripherals compatible with S Series ultrasound system (Continued)

Black & white printer cable	3.3 ft/1 m
Black & white printer USB cable	10.8 ft/3.3 m
Color printer	—
Color printer power cable	3.3 ft/1 m
Color printer video cable	6.0 ft/ 1.8 m
V-Universal Stand	—
Petite mouse	6 ft /1.8 m
Power cord (system)	10 ft/3 m
Power supply with DC cable	6.8 ft/2 m
Power supply AC cable	39 in/1 m
PowerPark	—
Triple Transducer Connect	—
USB keyboard	—
USB wireless adapter	—

### Guidance and manufacturer's declaration

**WARNING:** Other equipment, even equipment that complies with CISPR emission requirements, can interfere with the S Series ultrasound system.

The S Series ultrasound system uses an IEEE 802.11 receiver and transmitter that utilizes the ISM frequency band from 2.412 to 2.4835 GHz and implements three methods of transmission:

- IEEE 802.11b with Complementary Code Keying (CCK), Differential Quaternary Phase Shift Keying (DQPSK), and Differential Binary Phase Shift Keying (DBPSK) at 17 dBm
- IEEE 802.11g with Orthogonal Frequency Division Multiplexing (OFDM) at 14 dBm
- IEEE 802.11n with OFDM at 14 dBm

## Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The S Series ultrasound system is intended for use in the electromagnetic environment specified below. The customer or the user of the S Series ultrasound system should assure that it is used in such an environment.

<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment</b>
RF emissions CISPR 11	Group 1	The S Series ultrasound system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The S Series ultrasound system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The S Series ultrasound system is intended for use in the electromagnetic environment specified below. The customer or the user of the S Series ultrasound system should assure that it is used in such an environment.

<b>Immunity Test</b>	<b>IEC 60601 Test Level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment</b>
Electrostatic Discharge (ESD) IEC 61000-4-2	±6.0KV contact ±8.0KV air	±6.0KV contact ±8.0KV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient burst IEC 61000-4-4	±2KV for power supply lines ±1KV for input/output lines	±2KV for power supply lines ±1KV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity (Continued)

The S Series ultrasound system is intended for use in the electromagnetic environment specified below. The customer or the user of the S Series ultrasound system should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Surge IEC 61000-4-5	±1KV line(s) to line(s) ±2KV line(s) to earth	±1KV line(s) to line(s) ±2KV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5s	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the S Series ultrasound system requires continued operation during power mains interruptions, it is recommended that the S Series ultrasound system be powered from an uninterruptible power supply or a battery.
Power Frequency Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity (Continued)

The S Series ultrasound system is intended for use in the electromagnetic environment specified below. The customer or the user of the S Series ultrasound system should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the S Series ultrasound system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended Separation Distance</p> $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 2.3 \sqrt{P}$ <p>800 MHz to 2,5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (<math>W</math>) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p>

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity (Continued)

The S Series ultrasound system is intended for use in the electromagnetic environment specified below. The customer or the user of the S Series ultrasound system should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Radiated RF IEC 61000-4-3 (continued)			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following symbol: 

*Note:  $U_T$  is the AC mains voltage prior to application of the test level.*

*At 80 MHz and 800 MHz, the higher frequency range applies.*

*These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.*

- Field strengths from fixed transmitters such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SonoSite ultrasound system is used exceeds the applicable RF compliance level above, the SonoSite ultrasound system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the SonoSite ultrasound system.
- Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**FCC Caution:** Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:


- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

## Immunity testing requirements

The S Series ultrasound system complies with the essential performance requirements specified in IEC 60601-1-2 and IEC 60601-2-37. Results of immunity testing show that the S Series ultrasound system meets these requirements and is free from the following:

- Noise on a waveform or artifacts or distortion in an image or error of a displayed numerical value that cannot be attributed to a physiological effect and that may alter the diagnosis
- Display of incorrect numerical values associated with the diagnosis to be performed
- Display of incorrect safety related indications
- Production of unintended or excessive ultrasound output
- Production of unintended or excessive transducer assembly surface temperature
- Production of unintended or uncontrolled motion of transducer assemblies intended for intracorporeal use

## Labeling symbols

Symbol	Definition
	Temperature limitation
	Atmospheric pressure limitation
	Humidity limitation
	RF transmitter
	Prohibited action: Do not push
	Mandatory action: Transport position required
	Manufacturer

Authorized representative in European Community

Attention, see the user guide

Follow instructions for use.

## Standards

### Electrical safety standards

IEC 60601-1:1988, International Electrotechnical Commission, Medical Electrical Equipment-Part 1. General Requirements for Safety.

IEC 60601-1:2005, Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.

AAMI / ANSI ES 60601-1:2005, Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.

IEC 60601-1-1:2000, Medical Electrical Equipment-Part 1-1. General Requirements for Safety-Section 1-1. Collateral Standard. Safety Requirements for Medical Electrical Systems.

IEC 60601-2-37:2001, International Electrotechnical Commission, Particular Requirements for the Safety of Ultrasonic Medical Diagnostic and Monitoring Equipment.

IEC 60601-2-37:2007, Medical electrical equipment - Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment.

CAN/CSA C22.2, No. 60601-1, Canadian Standards Association, Medical Electrical Equipment-Part 1. General Requirements for Safety.

CAN/CSA C22.2, No. 60601-1:08, Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance.

IEC 61157, International Electrotechnical Commission, Standard Means for the Reporting of the Acoustic Output of Medical Diagnostic Ultrasonic Equipment.

### EMC standards classification

IEC 60601-1-2:2007, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests.

CISPR 11, International Electrotechnical Commission, International Special Committee on Radio Interference. Industrial, Scientific, and Medical (ISM) Equipment - Radio-Frequency Disturbance Characteristics-Limits and Methods of Measurement. Classification for the ultrasound system, docking system, accessories, and peripherals when configured together: Group 1, Class A.

## Acoustic standards

NEMA UD 2-2004, Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment.

NEMA UD 3-2004, Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment, American Institute of Ultrasound in Medicine.

## Biocompatibility standards

AAMI/ANSI/ISO 10993-1, Biological evaluation of medical devices - Part 1: Evaluation and testing (2009).

AAMI/ANSI/ISO 10993-5, Biological evaluation of medical devices - Part 5: Tests for In Vitro cytotoxicity (2009).

AAMI/ANSI/ISO 10993-10, Biological evaluation of medical devices - Part 10: Tests for irritation and delayed-type hypersensitivity (2002).

AAMI/ANSI/ISO 10993-11, Biological evaluation of medical devices - Part 11: Tests for systemic toxicity (2006).

AAMI/ANSI/ISO 10993-12, Biological evaluation of medical devices - Part 12: Sample preparation and reference materials (2007).

## Airborne equipment standards

RTCA DO-160E, Radio Technical Commission for Aeronautics, Environmental Conditions and Test Procedures for Airborne Equipment, Section 21.0 Emission of Radio Frequency Energy, Category B. 118.

## DICOM standard

NEMA PS 3.15, Digital Imaging and Communications in Medicine (DICOM)-Part 15: Security and System Management Profiles.

## HIPAA standard

Health Insurance and Portability and Accountability Act, Pub.L. No. 104-191.

45 CFR 160, General Administrative Requirements.

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