

---

FUJIFILM SonoSite, Inc.  
21919 30<sup>th</sup> Drive SE  
Bothell, WA 98021-3904

**For release to:**

- SonoSite Internal
  - Authorized Repair Facilities
  - Authorized Dealers/Distributors
  - SonoSite Customers
  - SonoSite.com
- 

**S-Series Ultrasound System (P/N: P09417-xx), X-Porte (P/N P17860-xx, P17870-xx)  
Compatible with Sanicloth AF3**

**Problem:**

Request for new disinfectant, Sanicloth AF3, to be tested for chemical compatibility on the S-Series and X-Porte Ultrasound System and transducers.

**Solution:**

Sanicloth AF3, has been tested and approved for the S-Series and X-Porte Ultrasound System and the following probes: C11x, C60x, HFL38x, ICTx, L25x, L38xi, and P10x.

This service technical bulletin serves as an interim communication tool for customers with an immediate need for this disinfectant while the website is being updated.

If there are any questions please contact FUJIFILM SonoSite Technical Support:

Phone: 877.657.8118 (U.S. and Canada)  
1.425.951.1330 (Worldwide)

Email: [service@sonosite.com](mailto:service@sonosite.com)