

User Guide Supplement #P09053-01

- *M-Turbo Ultrasound System User Guide*
- *S Series Ultrasound System User Guide*
- *MicroMaxx Ultrasound System User Guide*
- *TITAN Ultrasound System User Guide*
- *iLook Personal Imaging Tool User Guide*
- *SonoSite Ultrasound System User Guide (180PLUS and SonoHeart ELITE)*

Introduction

This document provides additional safety information regarding the following products:

- M-Turbo™ ultrasound system
- S Series™ ultrasound system
- MicroMaxx® ultrasound system
- TITAN® ultrasound system
- iLook® ultrasound system
- 180PLUS™ ultrasound system
- SonoHeart® ELITE ultrasound system

Safety

Hazardous materials

WARNING: | The liquid crystal display (LCD) contains mercury. Dispose of the LCD properly in accordance with local regulations.

Caution:

Federal (United States) law restricts this device to sale by or on the order of a physician.



SonoSite, Inc.

21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or
1-425-951-1200
F: 1-425-951-1201

SonoSite Ltd

Alexander House
40A Wilbury Way
Hitchin, Herts
SG4 0AP UK
T: +44-1462-444800
F: +44-1462-444801

03/2008
Copyright 2008 by
SonoSite, Inc.
All rights reserved.

P09053-01

